


What is χ ?

This MRI-based technology evaluates local magnetic field distortions to **characterize local tissue properties** in the liver. Each tissue has **intrinsic magnetic properties** that are different from one tissue to another. A tissue can be, to different extent, either diamagnetic, paramagnetic or ferromagnetic, leading to variations of local magnetic susceptibility.

For instance, **liver fibrosis** is characterized by the presence of Collagen type I & III in the extracellular matrix. This collagen is strongly diamagnetic leading to **negative χ** values while healthy liver tissues are slightly paramagnetic leading to positive QMS values. (figure 1)

Figure 1: Liver χ values after compensation for fat and iron confounding effects



SoQut
imaging



SoQut
imaging

ELiBio®
Digital Liver Biopsy

SoQut Imaging revolutionize diagnosis and management of Metabolic Associated Steatosis Hepatitis (MASH) through innovative, **non-invasive MRI** software-based technology.

Turn your **MRI** into a **biopsy**

www.SoQut-Imaging.com

Your **All-in-one software** solution for MALFD/MASH patient management

ELiBio[®] is an advanced MRI-based imaging software allowing **quantitative liver biomarkers** assessment such as Fat and iron contents. Most importantly, thanks to our unique quantitative **magnetic susceptibility** (χ - CHI) mapping technology, we can assess the presence and severity of **liver fibrosis**.



Compatible with most clinical MRI device (1.5T, 3T)



Quantitative metrics: Fat (%), Iron (mg/g), fibrosis (χ in ppm), visceral fat (%), lean mass (%) & subcutaneous fat (%)



3D mappings



Whole liver & spatial heterogeneity assessments



Smooth integration into radiological workflow: 20s, breath-holding modified mGRE T1 sequence



Time-to-result : < 10 mins




Clinical Case

- Men, 66 yrs, overweight, type II diabetes and controlled hypertension.
- Cytolysis since 1996, following HBV vaccination. Cytolysis peaks occurred several times, with progressive normalization each time

Clinical analysis

- Fibroscan[®]* LSM suggesting liver cirrhosis (14.4kPa-M probe, IQR/med 14%)
- MRI (without ELiBio) showed cirrhotic-type liver dysmorphism, with diffuse steatosis overload.
- The biopsy showed no fibrosis (SAF S2 A3 F0/1c)

ELiBio analysis



This patient Fibroscan and biopsy results are contradictory.

Using ELiBio, the patient has an average magnetic susceptibility representative of a healthy liver but with a strong spatial heterogeneity.

ELiBio Chi mapping highlights focal areas of fibrosis at the junction of the right & left lobes and in the anterior part of the left lobe.

Both Fibroscan and biopsy analysis were true and are confirmed by ELiBio analysis.